# Shoebox Spectroscope

The shoebox spectroscope takes the light from a slit at one end and passes it through a grating to be diffracted or be broken down into the colors it is made of (the spectrum).

**To make this simple instrument you will need:**

1. A shoebox or any thing similar to a shoe box with a long length and small width preferably made of cardboard.
2. A flashlight
3. An index card cut in to 2 parts
4. A diffraction grating slide
5. Scissors or a knife or any puncturing tool intended for cutting holes in the shoebox
6. Glass lens (optional) for clearer appearance of diffracted light
7. Tape or glue
8. A fairly dark environment with a nice reflecting wall

## PROCEDURE

1. Cut one hole, the size of the diffraction grating, on one end of your shoebox and tape the diffraction grating over the hole as shown.
2. Cut another hole the same size as the first on the opposing side of the shoebox
3. Create a cracklike opening (or slit) between the split pieces of the index card. Tape them in place making sure that the slit is centered on the end of the shoebox (see 3<sup>rd</sup> photograph).
4. Put the lid on the shoebox.
5. Turn on the flashlight and direct the light into the slit opening.
6. Go to a fairly dark environment with the flashlight and your homemade spectroscope.
8. Aim the diffraction grating window toward the wall and see the diffracted light projected on the wall.
9. Put a lens in front of the diffraction grating to focus the image

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